

T. Berg¹, F. Fuest¹, M. Schütte¹, J. Pareja², T. Li², A. Dreizler², B. Böhm²

¹LaVision GmbH, Göttingen, Germany ²TU Darmstadt, Germany
ffuest@lavisoin.de

Motivation

- ▶ Auto-ignition (AI) of turbulent jets with a microwave-plasma heated co-flow.
- ▶ Turbulent flow phenomena are three-dimensional in nature.
- ▶ Unsteady auto-ignition process requires high temporal resolution.
- ▶ **Tomographic OH-LIF** imaging as an approach to yield three-dimensional OH distributions.
- ▶ **High-speed** tomographic LIF reveals information necessary to capture and follow auto-ignition processes.

Microwave plasma heater burner

The burner is designed for the investigation of auto-ignition of fuel jets propagating into a hot co-flow at high velocity and turbulence.

- ▶ Co-flow: Air, 82 mm nozzle diameter
- ▶ Jet: CH₄, 6 mm nozzle diameter

Table 1: Operating conditions

$Re_{co-flow}$ (-)	$T_{co-flow}$ (°C)	$U_{co-flow}$ (m/s)	Re_{jet} (-)	$T_{co-flow}$ (°C)	U_{jet} (m/s)
10,000	1050	25	3200	520	51
10,000	1100	27	3700	530	60




Figure 1: Plasma heater co-flow burner.

Experimental setup




Figure 2: Schematic setup of the 8-view high-speed camera system, laser and burner.




Figure 3: Photography of the setup with burner in operation.

Detection

- ▶ Four intensified high-speed cameras
- ▶ 10 kHz frame rate
- ▶ UV lens f/8 with image doubler
- ▶ 8 simultaneous views

Excitation

- ▶ High-speed pump laser at 10 kHz
- ▶ Frequency-doubled output of a high-speed dye laser tuned to excite the Q1(6) transition ($\lambda = 283.01$ nm) of the $A^2\Sigma^+ \leftarrow X^2\Pi$ ($v' = 1, v'' = 0$) band of hydroxyl radicals
- ▶ 0.2 mJ at sample volume of $2.0 \times 1.5 \times 0.9$ cm³

Tomographic reconstruction

- ▶ **Simultaneous Multiplicative Algebraic Reconstruction Technique (SMART)**
- ▶ 100 iterations
- ▶ Volume smoothing after each iteration with $3 \times 3 \times 3$ voxel filter with smoothing factor of 0.5
- ▶ Computational time for 8 views @ 16 cores (3.10 GHz, 128 GB RAM):
 - ▶ 40 s for 32M voxel of $95^3 \mu\text{m}^3$, no binning
 - ▶ 4 s for 4M voxel of $190^3 \mu\text{m}^3$, 2×2 binning




Figure 4: Iso-surfaces of OH-LIF intensity. $Re_{jet} = 3200$, $T_{co-flow} = 1050$ °C




Figure 5: Reconstructed z-planes of OH-LIF from a flame kernel. $Re_{jet} = 3200$, $T_{co-flow} = 1050$ °C

Time-resolved 3D-imaging

- ▶ Time-series of reconstructed tomographic OH-LIF images for different operating conditions
- ▶ Re -number of the jet and co-flow temperature have a strong effect on the auto-ignition of turbulent CH₄ jets




Figure 6: Propagation of auto-igniting kernels: overlay image of multiple time steps. $Re_{jet} = 3700$, $T_{co-flow} = 1100$ °C




Figure 7: Time series of an auto-igniting CH₄ flame. $Re_{jet} = 3200$, $T_{co-flow} = 1050$ °C

Conclusions

- ▶ Applications of high-speed tomographic OH-LIF successfully proven on turbulent auto-igniting jet flames.
- ▶ Reconstructions of the case with $Re_{jet} = 3700$ and $T_{co-flow} = 1100$ °C showed the downstream-propagation and shape-evolution of auto-igniting flame kernels.
- ▶ 3D tomographic reconstructions of the case with $Re_{jet} = 3200$ and $T_{co-flow} = 1050$ °C showed spatial and temporal propagation of auto-igniting flame base.
- ▶ Higher framing rates of the tomographic system are desired for future investigations of the temporal evolution of auto-igniting flames and kernels.